### 

## : HEINE Binocular Loupes and Headlights

Perfect view for highest accuracy during all procedures. As individual as your daily work.

## **EXPERIENCE THE DIFFERENCE**



# : Experience HEINE<sup>®</sup> ML4 LED HeadLight


### Perfect fit. Perfect view. During long procedures.

Specially designed for healthcare professionals who need a perfect view during long surgical operations or treatments. The headband offers multiple adjustments points and soft padding ensuring highest comfort and a firm fit.



### **TOPSELLER IN ALL SURGICAL SPECIALTIES**

LED NOW IN HEINE QUALITY.





# : Experience HEINE® LoupeLight 2

Brighter. Lighter. Improved functionality.

# **MOST TRUSTED BY DENTISTS AND ENT SPECIALISTS**



> For ordering information, please see the last page or visit www.heine.com

# : HEINE® MicroLight 2

Ultralight. Universal. Without magnification.







# : Experience HEINE® Binocular Loupes

High resolution. Crisp image. Various magnifications.

Your optimal loupe is as individual as your work. HEINE's HIGH TECH optics offer a large variety of loupes to fulfill your individual needs (see page 6).

Experience the difference of getting less tired – even during long procedures through

- Crisp, high-resolution image
- Large field of view
- Homogenous magnification
- Lightweight
- Multi-coated, achromatic technology (HR) or
- Prismatic for high resolution (HRP)

Unique HEINE i-View Loupe Mount Adjustable viewing angle ensures a coaxial, shadow-free image.

Flip-up function. Ability to flip optics up for unmagnified views and overview.

## **OPTIMAL MAGNIFICATION**





Without light and without magnification



With a comparable competitor product (with illumination and with magnification)



With HEINE 2.5 x

**Binocular Loupes** 

and LoupeLight 2

LED NOW IN HEINE QUALITY.

# : Experience unique HEINE LEDHQ

Bright. Homogenous. True colour rendering.

LED is not just LED. Experience five important differences of HEINE LED<sup>HQ</sup> quality for perfect illumination during your daily work.



### **OPTIMUM BRIGHTNESS**

55,000 – 65,000 Lux, the proven optimum for our Loupelights and Headlights, will allow you to discern even slightest abnormalities.



### EDGE TO EDGE HOMOGENEITY

A coaxial, absolutely bright and uniform light spot.



### TRUE COLOUR RENDERING

Comparable with daylight, indicated by a Colour Rendering Index (CRI) of 90: Red is red and blue is blue.



### IDEAL TEMPERATURE MANAGEMENT

The compact, ergonomic design with a heat-conducting foil and an aluminium heat sink ensures optimal LED performance and life time.



### UNLIMITED WORKING LIFE

The Integrated HEINE LEDs<sup>HQ</sup> in HEINE Headand Loupelights have virtually unlimited working life and are maintenance-free.



## **PERFECT ILLUMINATION**



### ... AND BRIGHT AND HOMOGENOUS ILLUMINATION AND TRUE COLOUR RENDERING



Without light and without magnification

With a comparable competitor product (with illumination and with magnification)



With HEINE 2.5x Binocular Loupes and LoupeLight 2

LED NOW IN HEINE QUALITY.

# : Experience your perfect combination

In only 4 steps. Choose your perfect configuration.





# : Experience true individuality

Ordering information.

## **PERFECT VIEW**



### Our most trusted sets:

### LOUPES AND LIGHTS

ML4 LED HeadLight with HR/HRP Binocular Loupes with i-View loupe mount and S-GUARD on Headband Professional L,

### mPack and plug-in transformer, cleaning fluid

Binocular Loupe 2.5x/420mm working distance	J-008.31.441
P Binocular Loupe 3.5x/420mm working distance	J-008.31.442
P billocular Loupe 3.5 x7 42011111 Working distance	_

NEW! LED LoupeLight 2 with HR/HRP Binocular Loupes with i-View loupe mount on S-FRAME, detachable yellow filter, mPack mini, mPack mini Belt Clin, E4-I ISB plug-in transformer, protective lenses, 2 sterilisable levers, retaining cord, cleaning fluid in case

with HR Binocular Loupe 2.5x/420mm working distance	C-008.32.451
with HRP Binocular Loupe 3.5x/420mm working distance	C-008.32.453

### LIGHTS ONLY

ML4 LED HeadLight on Headband Professional L with Wall Transformer EN 50 UNPLUGGED, 2x mPack UNPLUGGED	J-008.31.416
NEW! MicroLight 2 on lightweight headband with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer,	
retaining cord and cleaning fluid, in case	J-008.31.277

### LOUPES ONLY

Binocular Loupes Sets with i-View on S-FRAME, Protective lenses, 2 sterilisable swivel levers, retaining cord, cleaning fluid,

cleaning cloth, carrying case

with HR Binocular Loupe 2.5x/420mm working distance	C-000.32.356
with HRP Binocular Loupe 4x/340mm working distance	C-000.32.431



Find your individual combination with all the various loupes options on www.heine.com/binocular\_loupes\_headlights



#### HEINE QUALITY MADE IN GERMANY

In 1946 our founder, Helmut A. Heine started HEINE Optotechnik with one goal in mind; to build the highest quality diagnostic instruments available in the world. It was his belief that the highest quality instruments will always best support the health care professional in making the earliest and most accurate diagnosis. For over 70 years, our customers have been able to rely on this promise, and on their HEINE instruments for precise diagnostic performance, unique functionality and long lasting durability.

Today, as a global leader in primary diagnostic instruments with over 500 employees we remain committed to our founding principles. We are still a 100% family owned and managed company and care deeply about our employees, our customers, and the quality of our products. We continue to develop and manufacture HEINE instruments exclusively at our facilities in Germany, where we combine decades-long experience and craftsmanship with the most modern manufacturing technologies.

### www.heine.com



#### GERMANY

644

HEINE Optotechnik GmbH & Co. KG Kientalstr. 7 82211 Herrsching Tel. +49 (0) 81 52-38 0 Fax +49 (0) 81 52-3 82 02 E-Mail: info@heine.com www.heine.com

#### NORTH AMERICA

HEINE USA LTD. 10 Innovation Way Dover, NH 03820 Tel. (603) 7 42-71 03 Fax (603) 7 42-72 17 Toll Free (800) 367-4872 E-Mail: service@heine-na.com

#### AUSTRALIA

HEINE AUSTRALIA PTY. LTD. Unit 9, 98 Old Pittwater Road PO Box 7218 Warringah Mall NSW 2100 Tel. +61 (0) 2-99 38 95 00 Fax +61 (0) 2-99 39 23 05 E-Mail: info@heine.com.au

#### SWITZERLAND

HEINE (Schweiz) AG Tobeläckerstr. 9 CH-8212 Neuhausen am Rheinfall Tel. +41 (0) 52-6 72 22 66 Fax +41 (0) 52-6 72 63 77 E-Mail: info@heineschweiz.ch

We reserve the right to change specification without notice. @ 2017 HEINE Optotechnik GmbH & Co. KG. All rights reserved.

Your distributor

07/17. A-000.00.257 e

# "You can only treat what you see!"

Dr. Erik Vetter (MSc Master of Science Implantology, Diplomate ICOI Implantology, Periodontology, Aesthetic Dentistry) testing HEINE HRP Loupes.

