

# UBM Plus

The UBM Plus probe has a guarded tip to prevent damage to both the patient and the probe. The unit is used for anterior segment imaging. This unit is an incredibly portable device because the 48MHz probe plugs directly into any Microsoft Windows laptop or desktop computer. The UBM features an all-in-one probe design to eliminate signal loss and provide the sharpest images possible.

## Features

- High-definition imaging – Detailed structure definition, including cornea, iris, ciliary body, zonules, crystalline lens, and intraocular lens, as well as pathologies
- State-of-the-art probe design – Sharper, more focused images due to the elimination of signal loss
- Unsurpassed data analysis – Contains tools for measuring sulcus-to-sulcus, anterior chamber depth, positioning of intraocular lenses, and filtration angle of the eye
- Easy-to-set-up defaults, smooth customer interface, and automatic EMR exports that decrease examination time, increase patient throughput, and increase profit
- Adjustable video loop
- Full upgradeable software
- EMR compatible
- DICOM ready
- Multiple measuring velocities for ACD
- PDF reporting
- Fully adjustable gain-control before, during, or after scan
- Data archive – Network or external hard drive

